

The instructions below are provided as a reference tool for caregivers who are already familiar with the use of an insulin pump and with insulin therapy in general. Not all screens are shown. For more detailed information on the operation of the t:slim X2 insulin pump, please refer to its user guide.



1

Tap **0 grams** to enter the carbs for the user's bolus.

Note: If this area reads "units," the carb feature is turned off in the active profile.



2

Enter desired value. Be sure "grams" is displayed above keypad for food boluses. Tap to continue.



3


Tap **Add BG** to enter your blood glucose.

Note: If certain advanced features are being used, the current Dexcom G6 CGM reading may auto-populate to the bolus calculator. Please see the User Guide for more information.




4

Enter desired value. Be sure "mmol/L" is displayed above keypad when entering blood glucose (BG) values.



5

If a BG is entered that is below the target, but above 3.9 mmol/L, the option to reduce the bolus amount will appear. Tap to accept that reduction; otherwise, tap .



6

Tap the calculated units value to manually adjust recommended dose. Tap to continue.



7


Verify the dose and tap to confirm.

Note: Calculations above are based on preset insulin-to-carb ratios and correction factors, which may be set in Personal Profiles.



8

Tap to deliver the food bolus immediately. The BOLUS INITIATED screen will appear to confirm delivery has started.



9

To cancel the undelivered portion of the bolus, tap "X" next to BOLUS on the Home screen, then tap to confirm canceled bolus.

Important Safety Information: Tandem insulin pumps are intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The t:slim X2 insulin pump is approved for individuals 6 years of age and older. For Important Safety Information, please visit tandemdiabetes.com/safetyinfo.

© 2022 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care and t:slim X2 are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. Inc. Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries. All third-party marks are the property of their respective owners. ML-1004690_B

