

From the desk of:



Thank you for helping me manage my type 1 diabetes at school.



Start here



This student uses the **t:slim X2™ insulin pump** from **Tandem Diabetes Care** and we thought it would be useful for you to have some additional information.



The pump uses **Basal-IQ® technology**, which can help this student spend less time worrying about lows by turning insulin delivery on and off as often as every five minutes (based on Dexcom G6 continuous glucose monitoring).

Check this out



Your student is able to do the following:

	Always	Sometimes	Never
Accurately count carbs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Identify symptoms when experiencing hyperglycemia (high blood sugar) or hypoglycemia (low blood sugar)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Independently navigate the t:slim X2 insulin pump screens to acknowledge and resolve alerts and alarms	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Independently navigate the t:slim X2 insulin pump screens to deliver a food and/or correction bolus	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Change an infusion set independently	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Comfortably discuss their insulin pump and/or diabetes status with classmates	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Other details about your student and their care at school:



Please visit tandemdiabetes.com/back-to-school for more information on how you can help any students living with type 1 diabetes who use the t:slim X2 insulin pump. You can even test the pump's interface using our t:simulator™ app.



Basal Rate: A slow, continuous delivery of insulin which keeps blood sugar levels stable between meals and during sleep. A basal rate is measured in units per hour.

Bolus: A quick dose of insulin that is delivered to cover food consumed or elevated blood sugar.

Personal Profile: A personalized group of settings that define the delivery of basal and bolus insulin within specific time segments throughout a 24-hour period.

Insulin-to-Carbohydrate Ratio: The number of grams of carbohydrate that one unit of insulin will cover.

Type 1 Diabetes: A condition in which beta cells in the pancreas are destroyed, preventing the body from producing insulin. People with type 1 diabetes must use insulin to treat this condition. Formerly known as Juvenile Onset Diabetes or Insulin-Dependent Diabetes Mellitus, IDDM.



Learn more at

TANDEMDIABETES.COM/BASALIQ

The information provided here is intended to be complementary to any existing 504 or official school-care plan. It does not replace a 504 or the official school-care plan from our healthcare provider.



Basal-IQ technology is not a substitute for active diabetes management and does not prevent hypoglycemia in all scenarios. Basal-IQ technology will not be able to predict glucose levels and suspend insulin delivery if your CGM is not working properly or is unable to communicate with your pump. Always pay attention to your symptoms and blood glucose levels and treat according to your healthcare provider's recommendations.



TANDEM
DIABETES CARE

(877) 801-6901
tandemdiabetes.com

t:simulator App
A free virtual pump demo



Important Safety Information: RX ONLY. The t:slim X2 pump and Basal-IQ technology are intended for single patient use for the management of diabetes mellitus in persons six years of age and greater. The t:slim X2 pump and Basal-IQ technology are indicated for use with NovoLog or Humalog U-100 insulin. **t:slim X2 insulin pump:** The t:slim X2 insulin pump with interoperable technology is an alternate controller enabled (ACE) pump that is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in people requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. **Basal-IQ technology:** Basal-IQ technology is intended for use with a compatible integrated continuous glucose monitor (iCGM, sold separately) and ACE pump to automatically suspend delivery of insulin based on iCGM readings and predicted glucose values. The bolus calculator is indicated for the management of diabetes by people with diabetes by calculating an insulin dose or carbohydrate intake based on user entered data. Users of the t:slim X2 pump and Basal-IQ technology must: be able and willing to use the insulin pump, CGM, and all other system components in accordance with their respective instructions for use; test blood glucose levels as recommended by their healthcare provider; demonstrate adequate carb-counting skills; maintain sufficient diabetes self-care skills; see healthcare provider(s) regularly; and have adequate vision and/or hearing to recognize all functions of the pump, including alerts, alarms, and reminders. The t:slim X2 pump, and the CGM transmitter and sensor must be removed before MRI, CT, or diathermy treatment. Visit tandemdiabetes.com/safetyinfo for additional important safety information.

© 2021 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, Basal-IQ, t:simulator, and t:slim X2 are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play and the Google Play logo are trademarks of Google LLC. All other third-party marks are the property of their respective owners. ML-1008360_A