



Alerts & Reminders

t:slim X2	Tandem Mobi
Options > My CGM	Settings > Alerts & Sounds
High Alert	<input type="checkbox"/> ON, alert me above: _____ <input type="checkbox"/> OFF
Low Alert	<input type="checkbox"/> ON, alert me below: _____ <input type="checkbox"/> OFF
Rise Alert	<input type="checkbox"/> ON <input type="checkbox"/> OFF
Fall Alert	<input type="checkbox"/> ON <input type="checkbox"/> OFF
Out of Range Alert	<input type="checkbox"/> ON, _____ minutes <input type="checkbox"/> OFF
Options > My CGM	Settings > CGM
Dexcom G6	Transmitter ID _____
Dexcom G7	Pairing Code _____
FreeStyle Libre 2 Plus	A new sensor is required when starting a t:slim X2 insulin pump. Not compatible with Tandem Mobi.

t:slim X2	Tandem Mobi
Options > My Pump	Settings > Alerts & Sounds
Low Insulin	<input type="checkbox"/> ON, _____ units <input type="checkbox"/> OFF
Auto-Off	<input type="checkbox"/> ON, _____ hours <input type="checkbox"/> OFF
Low BG	<input type="checkbox"/> ON, _____ mg/dL Remind me after _____ <input type="checkbox"/> OFF
High BG	<input type="checkbox"/> ON, _____ mg/dL Remind me after _____ <input type="checkbox"/> OFF
After Bolus BG	<input type="checkbox"/> ON, _____ mg/dL Remind me after _____ <input type="checkbox"/> OFF
Missed Meal Bolus	<input type="checkbox"/> ON _____ days <input type="checkbox"/> OFF

Control-IQ Technology

Load

t:slim X2	Tandem Mobi
Options > My Pump	Settings > Pump
Weight	
Total Daily Insulin	

t:slim X2	Tandem Mobi
Options > Load	Actions
Site Reminder	<input type="checkbox"/> ON, every ___ days <input type="checkbox"/> OFF
Fill Cannula	

Note: Use last 14-day average. Average can be found:
1. In the Overview Report of the Tandem Source platform (source.tandemdiabetes.com).
2. In your t:slim X2 insulin pump under Pump History

Note: In order to access the Fill Cannula amount, you will need to stop insulin delivery. DO NOT stop insulin delivery until you are ready to disconnect from your pump.

Important Safety Information: RX ONLY. Indications for Use: Tandem Mobi system: The Tandem Mobi insulin pump with interoperable technology (the pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription. The pump is indicated for use in individuals 6 years of age and greater. **t:slim X2 insulin pump:** The t:slim X2 insulin pump with interoperable technology is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in people requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient use. The pump is indicated for use with NovoLog or Humalog U-100 insulin. The pump is indicated for use in individuals 6 years of age and greater. **Control-IQ technology:** Control-IQ technology is intended for use with compatible integrated continuous glucose monitors (iCGM, sold separately) and alternate controller enabled (ACE) pumps to automatically increase, decrease, and suspend delivery of basal insulin based on iCGM readings and predicted glucose values. It can also deliver correction boluses when the glucose value is predicted to exceed a predefined threshold. Control-IQ technology is intended for the management of Type 1 diabetes mellitus in persons 6 years of age and greater. Control-IQ technology is intended for single patient use. Control-IQ technology is indicated for use with NovoLog or Humalog U-100 insulin.

WARNING: Control-IQ technology should not be used by anyone under the age of 6 years old. It should also not be used in patients who require less than 10 units of insulin per day or who weigh less than 55 pounds.

Control-IQ technology is not indicated for use in pregnant women, people on dialysis, or critically ill patients. Do not use Control-IQ technology if using hydroxyurea. Users of a Tandem insulin pump and Control-IQ technology must use the insulin pump, CGM, and all other system components in accordance with their respective instructions for use; test blood glucose levels as recommended by their healthcare provider; demonstrate adequate carb-counting skills; maintain sufficient diabetes self-care skills; see healthcare provider(s) regularly; and have adequate vision and/or hearing to recognize all functions of the pump, including alerts, alarms, and reminders. The Tandem pump and the CGM transmitter and sensor must be removed before MRI, CT, or diathermy treatment. Visit tandemdiabetes.com/safetyinfo for additional important safety information.

© 2024 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, the Tandem logo, Tandem Source, Control-IQ, t:slim X2, and Tandem Mobi are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. Dexcom, Dexcom G6, Dexcom G7, and any related logos and design marks are either registered trademarks or trademarks of Dexcom, Inc. in the United States and/or other countries. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott and used with permission. All third-party marks are the property of their respective owners. ML-1013483_B



(877) 801-6901
tandemdiabetes.com

FOR USE IN THE
UNITED STATES ONLY