






Insulin delivery starts once a Personal Profile is activated and can be stopped at any time. Use the following steps to stop all insulin delivery.



Note: These instructions are provided as a reference tool for pump users and caregivers who are already familiar with the use of an insulin pump and with insulin therapy in general. Not all screens are shown. For more detailed information on the operation of the Tandem Mobi system, please refer to its user guide.




Tap **Actions** from the Navigation bar.




Tap **Stop Insulin**.

Note: When insulin delivery is stopped, any active bolus or Temp Rate and Control-IQ technology adjustments are immediately stopped.



Use the smartphone's security feature to verify identity.

Note: The system uses device authentication to ensure data security and patient confidentiality.




Tap **Yes** to confirm.


Note: A red exclamation mark appears next to the insulin level indicator of the mobile app.

+


Control-IQ technology does not automatically resume insulin if it is stopped manually. Use the following steps to resume all insulin delivery, so Control-IQ technology can continue to automatically calculate adjustments.



Tap **Actions** from the Navigation bar.



Tap **Resume Insulin**.



Tap **Yes** to confirm.

🚩 **Note:** If authentication has expired, use the smartphone's security feature to verify identity.

🚩 **Note:** If a Temp Rate was active when insulin delivery was stopped, it will resume if time remains on the duration timer.

Temp Rate

A Temp Rate is used to temporarily increase or decrease the programmed basal rate (steady amount of insulin delivered throughout the day and night), which can be helpful for situations such as illness.


Because Control-IQ technology automatically adjusts basal delivery, it must be turned off before a Temp Rate can be started. **Flip to the back** to learn how to start a Temp Rate.



Temp Rate (continued)




Use the following step-by-step instructions to program the desired percentage and length of time for a temporary basal rate.



Tap **Actions** from the Navigation bar and then tap **Temp Rate** to continue. Use the smartphone's security feature to verify identity.




Tap **100%** to enter the desired percentage for the Temp Rate. The current programmed basal rate is 100%.




Using the on-screen picker, select desired percentage. Tap **Done**.


Note: The rate can be set from 0% to 250%. An increase is greater than 100% and decrease is less than 100%.



Tap **15 min** to program the desired length of time for the Temp Rate. Using the on-screen picker, select desired length of time. Tap **Done**.



Verify the Affected Time Segments and Adjusted Basal Rate(s) and tap **Start**.



Once enabled, an icon will appear next to the insulin level indicator and the basal rate delivery bar will be yellow.



To manually stop a Temp Rate, tap **Stop** on the Temp Rate screen and then tap **Yes**.

Note: If using Control-IQ technology, manually turn it back on after Temp Rate has stopped.



Need help? We have an extensive library of materials to help you with your pump.



Control-IQ technology does not prevent all highs and lows. You must still bolus for meals and actively manage your diabetes. Please visit tandemdiabetes.com/responsible-use for more info.

Important Safety Information: RX ONLY. Indications for Use: Tandem Mobi system: The Tandem Mobi insulin pump with interoperable technology (the pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription. The pump is indicated for use in individuals 6 years of age and greater. **Control-IQ technology:** Control-IQ technology is intended for use with compatible integrated continuous glucose monitors (iCGM, sold separately) and alternate controller enabled (ACE) pumps to automatically increase, decrease, and suspend delivery of basal insulin based on iCGM readings and predicted glucose values. It can also deliver correction boluses when the glucose value is predicted to exceed a predefined threshold. Control-IQ technology is intended for the management of Type 1 diabetes mellitus in persons 6 years of age and greater. Control-IQ technology is intended for single patient use. Control-IQ technology is indicated for use with NovoLog or Humalog U-100 insulin.

WARNING: Control-IQ technology should not be used by anyone under the age of 6 years old. It should also not be used in patients who require less than 10 units of insulin per day or who weigh less than 55 pounds.

Control-IQ technology is not indicated for use in pregnant women, people on dialysis, or critically ill patients. Do not use Control-IQ technology if using hydroxyurea. Users of a Tandem insulin pump and Control-IQ technology must use the insulin pump, CGM, and all other system components in accordance with their respective instructions for use; test blood glucose levels as recommended by their healthcare provider; demonstrate adequate carb-counting skills; maintain sufficient diabetes self-care skills; see healthcare provider(s) regularly; and have adequate vision and/or hearing to recognize all functions of the pump, including alerts, alarms, and reminders. The Tandem pump and the CGM transmitter and sensor must be removed before MRI, CT, or diathermy treatment. Visit tandemdiabetes.com/safetyinfo for additional important safety information.

© 2024 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, the Tandem logo, Tandem Mobi, and Control-IQ are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. All third-party marks are the property of their respective owners. ML-1013344_A



TANDEM
Diabetes Care

(877) 801-6901
tandemdiabetes.com

FOR USE IN THE
UNITED STATES ONLY

