

PATIENT NAME (FIRST, LAST):	TRAINER NAME (FIRST, LAST):	DATE (YYYY/MM/DD): ____/____/____
-----------------------------	-----------------------------	--------------------------------------

PERFORM THE FOLLOWING TASKS:

(Make sure you are disconnected from t:slim X2™ pump while performing these.)

1. Set the current **time** and **date**
2. **Fill the tubing** with a minimum of 10 units of insulin.
3. **Create & Name** _____ **a new Personal Profile**

Time	Basal	Correction Factor	Carb Ratio	Target BG

- Insulin Duration: _____ hours
- Carbohydrates: ON OFF
- Max Bolus: _____ units

4. **Activate** the new Profile
5. **Correction Bolus** – Enter BG of 9.0 mmol/L. **Deliver bolus.** (*cancel bolus*)
6. **Standard Food with Correction** –Bolus for 30 g Carb, BG of 11.7 mmol/L. **Deliver bolus.** (*cancel bolus*)
7. **Duplicate** your new Profile and call it “Travel”.
 Add a new time segment at _____ pm :
 Increase the Basal to _____ and decrease the Correction Factor to _____. **Do not Activate.**
8. **Standard Food Bolus** – Add 35 grams Carb - plus 4 grams Carb - plus 14 grams Carb.
 What is the total Carb bolus? _____ **Do not deliver.**

9. **Look up the last bolus.**
10. **Delete** the “Travel” Profile
11. **Stop** all insulin deliveries.
12. **Resume** all insulin deliveries.
13. **Change Cartridge, Fill Tubing, Fill Cannula**

OPTIONAL ACTIVITIES

- Extended Bolus** – Enter 40 g Carb.
Extend – 20% NOW, 80% LATER for 1 hr 30 min.
Deliver bolus. (*cancel the bolus*)
- Set Temp Basal Rate** – to 60% for 2 hours.
Stop the temp rate
- Deliver a Quick Bolus** of _____
 grams units