

NAME	TITLE	DATE (MONTH/DAY/YEAR)
------	-------	-----------------------







Please answer each question with your best response.

Section 1: Tandem t:slim X2 Insulin Pump

- What is the minimum amount of insulin needed in the cartridge AFTER the tubing has been filled to complete the load process?**
 - 25 Units
 - 50 Units
 - 75 Units
 - 95-120 Units
- What is considered a best practice for charging Tandem insulin pumps?**
 - Allow the pump battery to fully deplete and then charge to 100%.
 - Charge the pump for a short period of time (10 to 15 minutes) every day.
 - Charge the battery overnight.
 - Only charge the pump from a computer when disconnected from the infusion site
- True or False: After a cartridge change, the Insulin Level Indicator on the pump home screen displays an estimate of the amount of insulin in the cartridge as depicted by a + sign in front of a preset number. After ten (10) units have been delivered, the plus + sign will disappear, and the insulin level will display.**
 - True
 - False
- The fastest way to view the current active Personal Profile, the last bolus and if Basal-IQ or Control-IQ technology is turned ON is by tapping:**
 - Tandem Logo
 - Bolus icon
 - Battery level
 - Status, next to insulin level
- What will happen if the Auto-Off alarm is ON and set for 12 hours?**
 - All insulin delivery will be stopped if no interaction with the pump screen (i.e., tap 1-2-3 to unlock) after 12 hours
 - The pump and the mobile app will become unpaired if no interaction for 12 hours
 - Basal-IQ and Control-IQ are turned off if no CGM connectivity after 12 hours
- If a user taps “Bolus” but does not complete/deliver a bolus or does not use the back arrow to return to the home screen, what will happen?**
 - Basal delivery will stop
 - The user will get an Occlusion Alarm
 - The user will get an Incomplete Bolus Alert
 - The battery will drain quickly

7. Put the steps in order to properly complete the Cartridge Fill process:

(Illustrations do not include EVERY step)

Steps 1-6	Steps of Cartridge Fill
A	<p>While the needle is still in the cartridge, release the plunger and allow it to come to a neutral position.</p> 
B	<p>Fill the syringe with insulin and remove air bubbles from syringe</p> 
C	<p>Insert syringe into cartridge, maintain pressure on the plunger and fill cartridge slowly.</p> 
D	<p>After filling the syringe with insulin, insert needle (for the first time) in the white fill port of the cartridge.</p> 
E	<p>Remove syringe from cartridge, turn it with needle facing upward and press on plunger to slowly remove air</p> 
F	<p>Pull back plunger until it is fully retracted to remove air from cartridge</p> 

8. True or False: To cancel or stop a bolus, tap the red X to the right of “Bolus”. If bolus delivery has started, tap “yes” to confirm stopping the bolus. The “Units Requested” and the “Units Delivered” will be displayed.

- True
- False

9. A patient receives a second occlusion alarm shortly after the first, and all insulin delivery has stopped. What is the first action that should be taken?
- Ignore the alarms if glucose level is not >250 mg/dL
 - Check glucose level and change infusion set and cartridge as there may be a blockage in the system
 - Give an insulin injection immediately
 - Disconnect from pump and call HCP
10. The Tandem t:slim X2 is an updatable pump. What does this mean?
- Tandem will send a new insulin pump free of charge when requested by patient
 - The patient must call Tandem every 3 months to ask for software updates
 - Patients will be notified by email when new software is available with instructions on how to obtain it.
 - The pump keeps you updated on your current sensor glucose level
11. The Basal Limit is defaulted to ____ u/hr.
- 1u
 - 3u
 - 2u
 - 5u
12. A user can create up to ____ different Personal Profiles with up to ____ time segments in each profile.
- 15, 6
 - 3, 12
 - 6, 16
 - 2, 24

Section 2: CGM

13. Which is NOT an effective strategy to optimize connection between the pump and transmitter:
- Wear pump and transmitter on the same side of the body
 - Set Out of Range Alert to "On" for 20 minutes
 - Wear pump screen facing the body
 - Keep pump and sensor no more than 20 feet from each other
14. What does it mean if this icon is displayed on the pump or mobile app screen?



- A CGM sensor session is active, but transmitter and pump are out of range.
 - The CGM sensor session has expired
 - The CGM session is in startup period
 - A BG calibration must be entered
15. True or False: Tandem Automated Insulin Dosing (AID) algorithms will work without CGM data.
- True
 - False

Section 3: Basal-IQ Technology

16. Match the delivery statement on the left with the correct response on the right to make a correct statement about Basal-IQ technology functionality:

Resumes insulin delivery _____	A are indicated by lower half of diamond icon shaded red, black "S" in a red box beside the insulin level, and red bar on the CGM graph
Reduces the programmed basal rate _____	B Never
Stops basal delivery _____	C when either the actual sensor glucose value increases, or insulin has been suspended for 2 hours in any 2.5 hour period
Suspensions _____	D when sensor glucose is predicted to be <80 mg/dL in the next 30 minutes
During a Basal-IQ suspension, _____	E normal bolus deliveries continue, but extended boluses will stop.

17. When should a patient check their glucose with a fingerstick when using Basal-IQ technology?

- a. Every time the pump suspends insulin.
- b. Every time that they need to deliver a bolus.
- c. Only when the pump resumes insulin delivery.
- d. When the sensor reading does not match the clinical symptoms that the user is experiencing.

Section 4: Control-IQ Technology

18. Control-IQ technology will increase the basal rate:

- a. When the sensor glucose is predicted to be >160 mg/dL in the next 30 minutes.
- b. When the sensor glucose reaches 160 mg/dL.
- c. When sensor glucose is predicted to be >180 mg/dL in the next 30 minutes.
- d. When the sensor glucose reaches 180 mg/dL.

19. What is considered a "best practice" for using Sleep Activity?

- a. Use it for naps if they are at least 2 hours or longer
- b. Program a Sleep Schedule to automatically start every day.
- c. Use the Sleep Activity all day, every day.
- d. Manually start a Sleep Activity when you go to bed.


20. When Control-IQ technology is turned on, what is the target BG for all bolus calculations?

- a. 112.5 mg/dL
- b. 110 mg/dL
- c. 120 mg/dL
- d. 140 mg/dL

21. When Control-IQ technology is turned on, what basal rate will be delivered when glucose levels are predicted to be between 112.5-160 mg/dL?


- a. The Basal Limit will be delivered
- b. Programmed basal rate will always be delivered
- c. Programmed basal rate with a maximum of 3.0 u/hr will be delivered
- d. 50% of the programmed basal rate will be delivered

22. On the screen below, what does the circled icon in right corner next to the status icon indicate?



- a. The pump has stopped insulin delivery
- b. The user has initiated a correction bolus
- c. The insulin level is low
- d. Control-IQ is delivering an automatic correction bolus

23. On the picture below of the Control-IQ enabled pump bolus screen, why would the box on the right not have a glucose level displayed?



- a. The t:slim X2 insulin pump does not have a CGM signal
- b. There is no CGM trend arrow on the t:slim X2 insulin pump
- c. There is no CGM data being received by the pump
- d. All of the above

24. What is required for Control-IQ technology to work?

- a. Control-IQ technology turned On with Total Daily Insulin and Weight entered
- b. Complete Personal Profile with Carbs turned On
- c. Active CGM sensor session with a CGM value and trend arrow
- d. All of the above

25. In the picture below, why was basal delivery not stopped when the glucose level was below the red line?



- The sensor glucose levels were expected to rise
- The sensor was in the warmup period
- Control-IQ /Basal-IQ technology was not turned on, indicated by no grey diamond in the left corner
- The pump suspended, however the user resumed insulin

26. Match the Control-IQ Alert on the left with the correction notification screen on the right.

<div>CONTROL-IQ LOW ALERT</div>	<div>Control-IQ has increased your insulin, but your sensor readings remain above 200 mg/dL.</div> <div>Check your cartridge, tubing and site, and test your BG.</div> <div>CLOSE</div>
<div>CONTROL-IQ HIGH ALERT</div>	<div>Control-IQ has delivered the maximum allowable insulin in a 2-hour period.</div> <div>Make sure your Total Daily Insulin is correct in Control-IQ settings.</div> <div>CLOSE</div>
<div>MAX INSULIN ALERT</div>	<div>Control-IQ has predicted that you will drop below 70 mg/dL in the next 15 minutes.</div> <div>Eat carbs and test your BG.</div> <div>CLOSE</div>

27. Control-IQ technology should not be used by individuals who use less than ____ units of insulin per day or weigh less than ____ pounds (kg), which are the minimum inputs required to initiate Control-IQ technology and for it to operate safely.

- 5 units, 50 lbs. (22.7kg)
- 10 units, 100 lbs. (45.5kg)
- 10 units, 55 lbs. (25kg)
- 20 units, 75 lbs. (34kg)

28. What does the red shaded area on the pump home screen indicate?



- a. The basal rate was decreasing.
- b. The basal rate was increasing.
- c. The basal rate was stopped.
- d. An autocorrection bolus was delivered.

Section 5: Tandem t:connect™ mobile app

29. What must be done before the pump can be connected to the t:connect™ mobile app:

- a. Set up a t:connect account
- b. Enable Bluetooth in the t:slim X2 pump
- c. Enable Bluetooth on the mobile device
- d. All of the above

30. What is a best practice with the t:connect™ mobile app ?

- a. Force quit/close the app everyday
- b. Open the app only when a patient is at healthcare provider's office
- c. Interact with the app daily to keep data automatically uploading to the cloud
- d. Interact with the app every month to keep data automatically uploading to the cloud

31. Where is the first place a user should go for t:connect™ mobile app support?

- a. The t:slim X2 pump User Guide
- b. www.tandemdiabetes.com/mobilesupport
- c. Call Customer Technical Support
- d. The App Guide under Settings in the t:connect mobile application

Section 6: Tandem t:connect™ mobile app with Mobile Bolus

32. True or False? When using the mobile bolus feature, a bolus must be stopped/cancelled on the device where the bolus was initiated.

- a. True
- b. False

33. True or False? If the pump and mobile device lose Bluetooth connectivity AFTER a mobile bolus delivery has started, the bolus will automatically be canceled.

- a. True
- b. False

34. What security settings are necessary to be able to use the mobile bolus feature in the t:connect mobile app?

- a. Voice recognition technology
- b. Device passcode or facial identification
- c. Find my phone feature
- d. All of the above


35. To use mobile bolus feature (select all that are true):

- a. User must have a t:slim x2 with Control-IQ technology pump
- b. Pump and smartphone must be paired and in range
- c. Must have Wi-Fi or cellular connectivity to the app
- d. User must have a t:connect account
- e. Smartphone security feature must be enabled
- f. All of the above

36. The App Guide:

- a. Is found by tapping "Settings" then "Help," is always up to date, and is available offline
- b. Is a resource for app troubleshooting
- c. Provides steps for getting started, including Smartphone set up and app tour
- d. All of the above

37. True or False? The Correction Bolus toggle (circled below) allows the user to accept or decline a recommended correction bolus. If the glucose entered is <70 mg/dL, the Correction toggle will not be available.



- a. True
- b. False

NAME	TITLE	DATE (MONTH/DAY/YEAR)
------	-------	-----------------------

1. _____
2. _____
3. _____
4. _____
5. _____
6. _____
7. _____
8. _____
9. _____
10. _____
11. _____
12. _____
13. _____
14. _____
15. _____
16. _____
17. _____
18. _____
19. _____
20. _____
21. _____
22. _____
23. _____
24. _____
25. _____
26. _____
27. _____

28. _____
29. _____
30. _____
31. _____
32. _____
33. _____
34. _____
35. _____
36. _____
37. _____

Total of 37 Possible Points Passing Score is 37 (100%)	
_____ # correct	_____ %
ASSESSED BY NAME/TITLE _____	