


App Features	Tandem Mobi Mobile App for Tandem Mobi Pump	t:connect Mobile App for t:slim X2 Pump
Bolus From Your Smartphone	✓	✓
View Pump Data	✓	✓
Receive Alerts & Alarms [§]	✓	✓
Wirelessly Upload Pump Data	✓	✓
Stop & Resume Insulin Delivery	✓	 Only from Pump Screen
Enable Temp Basal Rates	✓	
Manage Personal Profiles	✓	
Manage CGM Sensor Integration	✓	
Program Sleep Schedules	✓	
Enable Exercise Activity	✓	
Extended Bolus Delivery	✓	

§ You cannot dismiss alerts and alarms from the t:slim X2 pump app and must do so on the pump itself, but once an alert or alarm has been dismissed on the pump it will also disappear from the mobile app.