


t:slim[®]

Insulin Pump



technical specifications



Rx Only

Physical Characteristics

Size

- 2.0" x 3.13" x 0.6"
(H x W x D)

Weight

- 3.95 ounces (with full 300 unit cartridge)

Color

- Black aluminum case with chrome trim

User Interface

- Bright, color touch screen

Screen

- Shatter-resistant glass
- 2.4" diagonal

Insulin

- Up to 300 units
- For use with Humalog[®] or Novolog[®] rapid-acting insulin

Battery

- Integrated rechargeable lithium polymer battery
- Typically lasts up to 7 days on full charge

Infusion Sets

- Standard luer-lock infusion sets

Safety Features

Insulin Delivery Mechanism

- Micro-Delivery[®] technology

Alarms

- Adjustable volume, visual, vibratory

Alerts/Reminders

- High BG (blood glucose), Low BG, After Bolus BG, Missed Meal Bolus
- Infusion Set Change

Watertight

- 3 feet for 30 minutes (IPX7 rating)

Insulin Delivery Settings

Personal Profiles

- 6 customizable profiles

Timed Settings

- 16 timed settings per personal profile

Adjustable Settings

- Basal rate, correction factor, carb ratio, target BG

Insulin Duration

- 2-8 hours in 1-minute increments

Bolus Insulin Delivery

Bolus Settings

- Grams of carbs or units of insulin

Bolus Options

- Standard Bolus, Correction Bolus, Extended Bolus, Quick-Bolus

Bolus Calculator

- Integrated calculator with numeric keypad

Bolus Increments

- 0.01 unit (at volumes greater than 0.05 units)

Insulin Sensitivity Factor

- 1:1 to 1:600 mg/dL (1 mg/dL increments)

Insulin-to-Carb Ratio

- 1:1 to 1:300 g (0.1 increments below 10 grams)

Maximum Bolus Size

- 25 units (with an option for an additional 25 units)

For a full list of technical specifications, please reference the t:slim® Insulin Pump User Guide in the customer support area of our website www.tandemdiabetes.com.

Basal Insulin Delivery

Basal Rate Increment

- 0.001 units at programmed rates ≥ 0.1 units/hr

Maximum Basal

- 15 units/hr

Temporary Basal Rate

- 15 minutes to 72 hours (range: 0%–250%)

Other

Pump History

- Up to 90 days of data (11,000 events)

Accessories

- Fun and fashionable pump accessories available – from classic leather designs to colorful protective cases

Warranty

- 4 years

© 2015 Tandem Diabetes Care, Inc. All rights reserved. t:slim, Micro-Delivery and Tandem Diabetes Care are registered trademarks of Tandem Diabetes Care, Inc. All other product names and trademarks are the property of their respective owners. U.S. and International Patents Pending. ML-000029_C.

In The Box Package Contents



t:slim Insulin Pump



t:slim Screen Protector (pre-applied)



t:clip™ Vertical Case



USB Cable



Wall Power USB Adapter



Car Power USB Adapter



t:slim User Guide



Cartridge Removal Tool